CENTRILAB

Laborzentrifuge
Lab Centrifuge

CENTRILAB

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The CENTRILAB is designed for the analysis of the behavior of suspensions. Target is the dimensioning of continuously working centrifuges.

The retention times in continuously working centrifuges are short. Simulation was difficult in the past because acceleration- and deceleration-times of lab-centrifuges were unknown and long. The real retention-time deviates from the preset retention-time because the product is exposed to a transient centrifugal-force during the acceleration which is not negligible.

Hence the CENTRILAB works with minimized and defined acceleration- and deceleration-times. It is possible to acceleration within 1 second to the maximum speed. To reduce the acceleration-time the acceleration has to be increased. But a high acceleration lets the suspension spill over and churn up the sedimentated solids. These effects will also distort results of tests.

To avoid spill over and churn up the CENTRILAB works with an especially developed soft-start-function. The soft-start-function is realized by a high performance drive system with a frequency converter. The acceleration time and retention time can be set on a touchscreen in steps of a tenth of a second.

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Screening Machines Process Equipment

circular and elliptical motion screens
double counterweight screens
multideck horizontal screen
round screens
jigs

Sampling Systems, Airtube Systems, Size Reduction Machines, Laboratory Equipment, Control Screening Machines and Automation

individual units and complete installations
for sample taking and preparation
airtube systems
jaw crushers
roller mills
hammer and hammer impact mills
eccentric vibrating mills and ball mills
control screening machines
analytical screening machines
splitter
testing drums
automation

Centrifuges

scroll-screen centrifuges
pusher centrifuges
sliding discharge centrifuges
vibratory centrifuges
decanter centrifuges